

RAPIDES

Alcohol-based rapid disinfectant for surfaces

- **Effective against a broad range of microorganisms and viruses: being able to destroy Gram-positive and Gram-negative bacteria (incl. *Mycobacterium terrae*) and fungi, deactivates viruses (HBV, HIV, HCV, Vaccinia, Herpes simplex, Adeno, Rota)**
- **Treated surfaces maintains bactericidal activity during 24 hours.**
- **Leaves no residue, it dries quickly.**

Areas of application

Is used for quick disinfection of alcohol-resistant surfaces. This product is used in medical institutions, laboratories, the food industry, in child care establishments and in all other areas with high hygienic requirements. Very good for rapid disinfection of medical equipment surfaces.

RAPIDES – is ready-to-use colorless liquid with alcohol-like smell

Active ingredients

Ethanol – 65,0%, Propan-2-ol – 4,7%,
Polyhexamethylene Biguanide Hydrochloride – 0,1%.

Packaging units

1 L - spray bottle; 5 L canister; wet-wipes 100 ut. or 150 ut.

Storage life

3 years in original package

Labelling according Regulation (EC) No 1272/2008 [CLP]

Signal word - Warning



Hazard statement(s)
H226 Flammable liquid and vapor.

Precautionary statement(s)
P210 Keep away from heat/sparks/open flames/hot surfaces. – No smoking.



Product is tested in accordance with

EN 1275:2006, EN 1276:2009, EN 13624:2013, EN 13697:2015,
EN 13704:2002, EN 13727:2012+A2:2015, EN 14348:2005,
EN 14349:2012, EN 14561:2006, EN 14562:2006, EN 14563:2009,
EN 14476:2013+A1:2015, ISO 22196:2011.

Efficiency

Microorganism	Contact time
Bacteria (incl. MRSA, <i>Pseudomonas aeruginosa</i>) and yeast (<i>C. albicans</i>)	30 sec
Dermatophyte- Trichophyton	1 min
Mycobacterium M.Tuberculosis	1 min
Mycobacterium M.Terrae	5 min
Enveloped viruses (incl. HIV, HBV, HCV, Vaccina, Herpes, Influenza)	30 sec
Unenveloped viruses (incl. HAV, Adeno и Polyo)	5 min
Mould A.niger	5 min

Recommended usage

Solution - Apply on the surfaces and after the time of effect wipe if necessary. The surfaces must be completely covered with the substance. Contact time is 30 – 60 seconds.

Wipes - Open the cover and pull out a wipe up to the perforation. To tear off the wipe, support the cover with the other hand. After tearing the wipe off, close the cap immediately to avoid drying of the wipes. Wipe the surfaces, make sure that the whole area is covered with the disinfectant. Contact time is 30 seconds.

Operating norm - 50 ml / m2.



Baltiachemi OÜ
Maagi 5
74114 Maardu, Estonia
T. +372 621 4694
info@baltiachemi.ee
www.baltiachemi.ee

